

National Science Foundation Plant Genome

Cereal Plant Transformation Workshop

Albert Kausch

University of Rhode Island

Microprojectile Bombardment

Kimberly Nelson-Vasilchik, Joel Hague, and Albert Kausch
University of Rhode Island



Microprojectile Bombardment

Kimberly Nelson-Vasilchik, Joel Hague, and Albert Kausch
University of Rhode Island

Microprojectile bombardment was developed by John Sanford and Ted Klein in the mid-1980's. Bombardment then proved as the protocol for the first achieved transformation of maize and other monocots resulting in fertile transformants. It was subsequently replaced by enhanced *Agrobacterium*-mediated protocols. However, this method is still important for several reasons and applications. First microprojectile bombardment, serves as a delivery method for transient expression assays. Second, this protocol has become increasingly important to the development of genome editing delivery.

Microprojectile Bombardment



**Biolistic PDS-1000/ He
Particle Delivery System
Catalogue number 165-2257
BioRad**

Biolytic PDS-1000/ He Particle Delivery System

Catalogue number 165-2257

BioRad

Macrocarrier Cover Lid

Microcarrier Launch Assembly

Bombardment Chamber Door with gasket

Rupture disc retaining cap

Brass Assembly Nest

Target Plate Shelf



Biologic PDS-1000/ He Particle Delivery System

3-way Helium Metering (Solenoid) Valve

Power cord

2.5 ft PEEK Tubing
to top of
bombardment
device

6 ft PEEK Tubing to
Helium tank


Reinforced Vacuum
Tubing to Vacuum Pump



Microprojectile Bombardment Gun Components



Gun Preparation



Autoclave removable metal gun parts (as shown, right), stopping screens, macrocarrier holders and forceps. Sterilize macrocarriers (henceforth called Flyers) in 70% ETOH for 1 hour, and then dry in hood. *Ensure these and all gun parts are completely dry.*

Gun Preparation



Surface Sterilize the laminar flow hood working surfaces and the gene gun (inside and out) with 70% ETOH.

Microprojectile Bombardment Gun Components



Macro carriers
Of "Flyer"



Rupture Discs

These are purchased with various psi rupture points. The PSI rupture point is related to the acceleration of the particle

Microprojectile Bombardment Gun Components



Various Gun Tools for making adjustments when required.

Gun Preparation

Surface sterilization of Flyers



Macrocarriers (also called “Flyers”)

Surface sterilization of Flyers



Sterilize macrocarriers (henceforth called flyers) in 70% ETOH for 1 hour

Surface sterilization of Flyers



Start surface sterilization in 70% ETOH for 1 hour.

Surface sterilization of Flyers



After surface sterilization in 70% ETOH for 1 hour, use sterile forceps to move the flyers to a sterile Petri dish placed toward the back of the hood to dry. Note: because they were sterilized in 70% ETOH the flyers may take 20-30 minutes to dry. They may be left in the biosafety cabinet until needed.

DNA Precipitation

Reagents: Prepared prior to procedure

2.5 M $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ (100ml):

Dissolve 36.76g in ~50ml ultrapure H_2O

Bring to a final volume of 100 ml

Filter sterilize or autoclave solution



DNA Precipitation

Reagents: Prepared prior to procedure



0.1 M Spermidine (1ml):

Mix 15.7 μ l spermidine with 984.3 μ l ultra pure water

Filter sterilize *Note: Spermidine when properly stored (2-8°C) is a solid. Melt spermidine in 65°C water. Spermidine will oxidize and turn yellowish and should therefore be stored under N_2 . In addition, keep bottle and a desiccant pack in a 50 ml Falcon tube. Mix a fresh 1 ml working solution each bombing campaign.*

DNA Precipitation

Particle Stock: Prepared prior to procedure



Prepare a gold stock solution. Transfer 60 mg of 0.6 μ gold to a 1.5 ml tube with 1 ml of 100% ETOH.

Note: 100% ethanol will hydrate quickly becoming 90-95%. So, Splurge and open a new bottle of ETOH. Using old ETOH is likely to cause clumping of gold particles and decrease the efficiency of the transformation.

DNA Precipitation

Prepare working solutions of gold particles and DNA. To maintain uniformity in bombardment parameters, only 3 bombardments should be performed per solution. Also, keep experiments under six to eight tubes



DNA Precipitation




Vortex the gold stock (60 mg/ml) for 1 min. to break up any gold aggregates.

DNA Precipitation



Sonicate the gold for 15-20 seconds. *Caution: Although it is important to find the “active spot” in the sonicator, leaving the gold in the sonicator for more than 20 seconds will overheat the gold particles*

DNA Precipitation




Invert the tube between sampling to prevent settling of gold particles, withdraw 35 μ l of suspended gold from the middle of the 1.5 ml Eppendorf tube and dispense it to the bottom of a sterile 1.5 ml tube. Invert the tube between each sample to ensure a uniform sample concentration of the gold, remember gold particles will settle quickly

DNA Precipitation




Using a centrifuge with a swinging bucket rotor, pellet the gold for 5 minutes at 500 RPM. (*This is a good time to prepare the spermidine.*)

DNA Precipitation




This is a good time to prepare the spermidine. Take an aliquot of meted spermidine from the stock and filter sterilize it into a sterile microfuge tube.

DNA Precipitation



Remove the ETOH and without disturbing the pellet (left), then carefully add 1 ml of sterile water without disturbing the pellet. If pellet is not disturbed proceed without an intermittent centrifugation step. Otherwise pellet the gold for 5 minutes at 500 RPM.

DNA Precipitation




Carefully, without disturbing the pellet, remove the H₂O (as close to close to dryness as pipeting will allow).

DNA Precipitation




Resuspend the microprojectiles in $9\mu\text{l}$ of DNA solution [$1\mu\text{g}/\mu\text{l}$]. Pipette up and down to facilitate re-suspension. *The microprojectiles should go into solution readily and not be clumpy. Multiple constructs can be added. Co-transformation can occur up to 70% of the time.*

DNA Precipitation




Add 220 μ l of sterile H₂O, and mix well.

DNA Precipitation



Add 250 μl of sterile CaCl_2 , and pipette up and down to mix well. *Calcium chloride salts the DNA out of solution.*

DNA Precipitation



Add 50 μ l of sterile spermidine, mix well and place on ice. *Spermidine tends to stabilize the DNA.*

DNA Precipitation



Using a Genie 2 vibrating mixing pad (left) mix this solution in refrigerator (4C) for 10 minutes. *A small amount of foam on the top of the solution indicates a good rate of vortexing. Note: This is a good time to start the flyers drying and start soaking the rupture discs*


Surface sterilization of Rupture discs



Note: This is a good time to start surface sterilization of the rupture discs in 70% ETOH. Place rupture discs in 70% ETOH for at least 20 minutes. The rupture discs should remain in 70% ETOH until used. They will be inserted into the gun rupture disc holder still wet with ETOH.

Soak rupture disks (1100psi) in 70% ETOH for at least 20 minutes prior to use. *Oversozaking may affect the integrity of the mylar material and rupture setpoint*

DNA Precipitation



Centrifuge in a swinging bucket rotor for 5 minutes at 500 RPM in a refrigerated centrifuge (*if available, this may help to reduce DNase activity.*)

DNA Precipitation



Remove the supernatant from the DNA/gold precipitate (as much as possible) and carefully add 600 μ l of 100% ETOH. Resuspend the pellet.. It may be necessary to break the pellet apart with a pipet tip.

DNA Precipitation



Remove the supernatant from the DNA/gold precipitate (as much as possible) and carefully add 600 μ l of 100% ETOH to each of the tubes. Recall that a freshly opened or properly stored bottle of 100% ETOH should be used. . Resuspend the pellet.. It may be necessary to break the pellet apart with a pipet tip.

DNA Precipitation



Spin for 5 minutes at 500 RPM in a refrigerated centrifuge, remove ETOH and re-suspend pellet in 38 μ l of 100% ETOH. The precipitate should readily go into solution. The particles are now ready for bombardment

Microprojectile Bombardment



Open helium tank supply valve, and adjust line pressure to ~1200 psi (slightly higher than rupture disc setpoint.)

Microprojectile Bombardment



Turn on the vacuum pump

Microprojectile Bombardment



Place flyers in holders. The flyer must be insert underneath the ridge inside the holder. This may take some patience. *Make several of these at a time. (Multiples of 3-number of bombs per tube.)*

Microprojectile Bombardment



Resuspend DNA/gold by sonicating for 10 seconds

Microprojectile Bombardment



Using a micropipeter, place 10 μl of the sonicated particle solution on center of flyer.
Allow the flyer to dry thoroughly.


Microprojectile Bombardment



Place stopping screen in the microcarrier launch assembly.

DO NOT FORGET TO THIS !!!!

Microprojectile Bombardment



Invert the flyer holder and insert into microcarrier launch body. Then screw flyer holder cap screw onto microcarrier launch body. *Note: the flyer holder can be inserted in different positions pending the placement of spacers. Start with the middle location unless a different location has been optimized for.*

Microprojectile Bombardment



Screw flyer holder cap screw onto micrcarrier launch body. *Note: the flyer holder can be inserted in different positions pending the placement of spacers. Start with the middle location unless a different location has been optimized.*

Microprojectile Bombardment




Unscrew rupture disk retainer cap from inside of the gun, and insert a wet rupture disk which has been soaking in 70% ETOH until use. Rupture discs are available with different rupture setpoints, and can be used in pairs to achieve desired setpoints.

Microprojectile Bombardment



Hand screw retainer cap back on and then tighten firmly (but hand tight...do not over tighten!) .

Microprojectile Bombardment




Inset loaded microcarrier launch body into the first slot from the top with the hole to the front (this distance may and is determined experimentally for each system)

Microprojectile Bombardment



Insert Petri dish platform into the third slot from the top

Microprojectile Bombardment




Insert Petri dish containing target tissue onto the platform

Microprojectile Bombardment



Close chamber door. Operator should wear safety goggles.

Microprojectile Bombardment



Begin drawing a vacuum by pressing upward on the vacuum button . Watch as the gauge goes up until it indicates 27-28”Hg.

Microprojectile Bombardment



Begin These valves control the Vacuum Flow Rate (left) and the vent Flow Rate (right)

Microprojectile Bombardment



Rupture disc will burst at 1100 psi Hg


Press down the button rapidly to the ‘Hold’ position . Press up on the ‘Fire’ position and monitor the gauge on top of the device. When the pressure exceeds rupture disk setpoint (in this case 1100 psi) the gun should fire. Note the actual fire pressure on the gauge located on the top of the gun
Operator should wear safety goggles.

Microprojectile Bombardment




After gun fires, release the fire button and vent the chamber. When the vacuum gauge read zero the chamber door may be opened and the sample removed and covered.

Microprojectile Bombardment



Remove microcarrier launch body and place on a sterile surface. Unscrew cap and remove flyer holder. Unscrew rupture disk retainer cap and remove spent rupture disk. Return to step 4 for additional transformation attempts.

Microprojectile Bombardment



Remove microcarrier launch body and place on a sterile surface. Unscrew cap and remove flyer holder. Unscrew rupture disk retainer cap and remove spent rupture disk. Return to step 4 for additional transformation attempts.

Shutting down the Gun Apparatus

Lastly, after all the bombardments have been completed:

1. Shut off the valve to the Helium tank.
2. Draw a vacuum on the gun chamber by closing the chamber door, and pressing the vacuum button.
3. Next press the vent button to deplete the helium from the line. This procedure should be repeated until the pressure reaches zero.
4. Vent the chamber and open the door.
5. Turn off the vacuum pump last (*Caution : DO NOT leave the chamber under vacuum and then turn off the pump as the pump oil will be sucked into the chamber*)

Thank You