Developing an Optical Microscopy Method for Identifying Microplastics in

Environmental Samples

UNIVERSITY

THE

Rowan Bixler¹, Mohammadjavad Hajirezaei², and Dr. Ryan Poling-Skutvik²

¹Department of Physics, Amherst College, ²Department of Chemical Engineering, University of Rhode Island




BACKGROUND

- Microplastic particles often result from the degradation of larger plastic materials [1].
- Their presence in coastal environments poses a growing threat to ecosystems and human health [1].
- Current detection methods are limited by lengthy processing times, high costs, low accessibility, and weak signal detection in complex samples [2].
- Developing an identification method using particle dynamics through microscopy imaging would provide valuable information to help track microplastics, identify their sources, assess exposure risks, and inform mitigation strategies [3].

How can particle dynamics and light scattering patterns observed with microscopy be used to identify the size and composition of microplastics?





 $\left|\Delta \hat{I}(q,t;\Delta t)\right|$


Image structure function


Time step, Δt (s)


RESULTS AND DISCUSSION

Microscope Imaging

§ Poly(methyl methacrylate) (PMMA) and polystyrene (PS) spheres (≈ 1 μm in diameter) were dispersed in water, and videos were recorded.









MATLAB Analysis

§ More than 60,000 video frames were analyzed.

PS particles in water





- The diffusivity increases from the bottom of the glass slide up to the middle of the sample, reaching a plateau where particle motion becomes
- § Amplitude Function A(q) and Scattering halo of an individual particle:




completely random.

§ Investigating the relationship between amplitude A(q) and the effective length scale of the particles








PS at different heights




 L_{eff} = 1.33









CONCLUSION & FUTURE WORK

- Microscopy imaging combined with DDM provides valuable insights into the physical properties of microplastics.
- DDM is a powerful technique for probing the dynamics of microplastics in water.

$$L_{eff} = \frac{1}{q_{pea}}$$

- Additional measurements with different plastic types are needed to broaden the analysis and validate trends.
- Measuring at multiple heights (z-positions) can improve the fitting.
- Analyze mixtures of various plastics

ACKNOWLEDGMENTS

of Rhode Island, and the Donald R. Wilson Jr. Family Foundation. Partial support for professional development activities was provided by EPSCoR Cooperative Agreement #OIA-2433276 and the RI Commerce Corporation through the Science and Technology Advisory Committee. Any opinions, findings, conclusions, or recommendations expressed in this material are those of the author(s) and do not necessarily reflect the views of the funding partners

REFERENCES


[1] Hossain, R,; Ghose, A,; Sahajwalla, V. Microplastics pollution in Australia: Sources, impacts, and mitigation strategies. In Case Studies in Chemical and Environmental Engineering, Volume 11. June 2025, 101036

[2] Ivleva, N. Chemical Analysis of Microplastics and Nanoplastics: Challenges, Advanced Methods, and Perspectives. Chemical

[4] Bayles, A. V.; Squires, T. M.; Helgeson, M. V. Dark-field differential dynamic microscopy. *Soft Matter*.2016, 12, 2440. [5] Giavazzi, F.; Brogioli, D.; Trappe, V.; Bellini, T.; Cerbino, R. Scattering information obtained by optical microscopy: Differential

DIFFERENTIAL DYNAMIC MICROSCOPY (DDM)

Differential Dynamic Microscopy (DDM) is an image intensity fluctuations [4,5].



 $D(q, \Delta t) \equiv \left\langle \left| \Delta \hat{I}(q, \Delta t) \right|^2 \right\rangle$ Image Structure Function

 The dynamic image structure function decomposes analytically into three sub-functions:

$$D(q, \Delta t) = A(q)[1 - g(q, \Delta t)] + B(q)$$

- where A(q) is Amplitude function and B(q) is the background noise.
- For Brownian motion, the Intermediate Scattering function (g):

 $g(q, \Delta t) = exp(-\Gamma(q)\Delta t) \longrightarrow \Gamma(q) = q^2 D_{self}$

How can we use DDM analysis to extract key information about different microplastic chemistries?

Diffusivity = $0.40 \pm 0.01 \,\mu\text{m}^2/\text{s}$